

CorrectJitter Manual

22/03/2022

sebastien.schaub@imev-mer.fr

I. Goal


Correct unstable illumination in transmission microscopy.

II. How to Use

1. Open an XYt stack.
2. Run the macro either:

- from MicaJ :


Click on "MICA" > "Correct Jitter"



- from Standard mode :

Click on "Plugins" > "Macros" > "Run..." and select the "CorrectJitter.ijm" in your computer.

3. If no ROI has been drawn, the macro shows a minimum projection and wait from you drawing a rectangle ROI



Institut de la Mer de Villefranche/Mer

181, Chemin du Lazaret,


06234 Villefranche-sur-Mer Cedex, France

Sebastien.Schaub@imev-mer.fr



Plateforme d'Imagerie par Microscopie (PIM)

- The macro creates a new XYt stack with homogeneous illumination.



Institut de la Mer de Villefranche/Mer

181, Chemin du Lazaret,

06234 Villefranche-sur-Mer Cedex, France



Sebastien.Schaub@imev-mer.fr